

PMP2-peripheral myelin protein 2 Gene

Reference: PTXBC034997

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FABP8; M-FABP; MP2; myelin P2 protein; peripheral myelin protein 2		
For research use only.			

GENE DESCRIPTION	peripheral myelin protein 2
GENE ID	5375
ORIGIN SPECIES	Human
ACCESSIONS	BC034997
VECTOR	pDONR223
NCBI SYMBOL	PMP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaacaaattctgggcacctggaaacttctcttagtgagaactttgacgattacatgaaagctctgggtgtgggttagccaccagaaaactggaaatttggccaaaccactgtgatcatcagcaagaaaggagatattataactatacgaactgaaagtaccttataaaatag
aaatctccttcaagctaggccagggaatttgaagaaccacagctgacaatagaagaccagagcatcgttaacctgcagagaggatcactgaatcaagtgcagagatgggatggcaagagacaaccataaagagaaagctagtgaatgggaaatggtagcggatgtaaaatgaa
gggcgtgggtgtgcaccagaatctatgagaaggctga