



Reference: PTXBC035007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM71A; GARI-L4; Golgi-associated Rab2B interactor-like 4; testis secretory sperm-binding protein Li 209a; family with sequence similarity 71 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 71, member A
GENE ID	149647
ORIGIN SPECIES	Human
ACCESSIONS	BC035007
VECTOR	pDONR223
NCBI SYMBOL	FAM71A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr