

FUNDC1-FUN14 domain containing 1 Gene

Reference: PTXBC035015

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FUN14 domain-containing protein 1; FUN14 domain containing 1		
For research use only.			

GENE DESCRIPTION	FUN14 domain containing 1
GENE ID	139341
ORIGIN SPECIES	Human
ACCESSIONS	BC035015
VECTOR	pDONR223
NCBI SYMBOL	FUNDC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgacccggaacccccctcccaagactatgaagtgatgacgactcttatgaagtggttgatttaactgagatgcaagaagacaccagtggtggaatcgagtggttgccacagttcgggacattgtagaaaaatactcagtagctaccagattgtaatgggtggcgttactggc
tgggtgagcaggtttctgtccagaagttggaaaactgcagcaactgcagtaggtgggtggtcttctcttcagattgctagtcatagtggctatgtgcagattgactggaagagagttgaaaaagatgtaataaagcaaaaagacagattaagaacgagcgaacaaagcagcacctg
aatcaacaatttaattgaagaagcaacagaattatcaagcagaacattgtgatatccagtggtgattgtggaggcttttgcctcgactgcactctaa