

RAB28-RAB28, member RAS oncogene family Gene

Reference: PTXBC035054

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB28, member RAS oncogene family; CORD18; ras-related protein Rab-28		
For research use only.			

GENE DESCRIPTION	RAB28, member RAS oncogene family
GENE ID	9364
ORIGIN SPECIES	Human
ACCESSIONS	BC035054
VECTOR	pDONR223
NCBI SYMBOL	RAB28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggactctgaggaggagagtcaggaccgcaactgaaatcgtgctgggggacggcgctccgggaagacctcttaactacgtgttttgcctcaagaactttgggaaacagtacaacaaactataggactggattctttttgagaaggataacattgccaggaaactgaat
gttacccttcaaatgggatagggaggcagacaataggaggcaaaatgttgataaataatctatggagcacagggagtcctcttggtatatgatattacaattatcaagctttgagaatttagaagattggtatactgtggtgaagaaagtgcagaggagtcagaaactcagccact
gggtgccttgtaggcaataaaattgattggagcatatgcgaacataaaacctgaaaaacacttacgggttttgcaggaataatggttttagtagccactttgtctcagccaagacaggagactctgtctctgtcttcagaaagtgtgctgctgaaatccttgggatcaaatatacaagca
gaaatagaacagtcacagaggggtggaaggcagatattgtaactacaaccaggaaacctatgtcaaggactgttaaccctcctagaagctctatgtgtgcagttcagtgga