

MXI1-MAX interactor 1 Gene

Reference: PTXBC035128

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAD2; MXD2; MXI; bHLHc11; max-interacting protein 1; MAX dimerization protein 2; Max-related transcription factor; class C basic helix-loop-helix protein 11; MAX interactor 1, dimerization protein		

For research use only.

GENE DESCRIPTION	MAX interactor 1
GENE ID	4601
ORIGIN SPECIES	Human
ACCESSIONS	BC035128
VECTOR	pDONR223
NCBI SYMBOL	MXI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagcacaatcaacttttctgcagaacgtgcagattctgtcaggccgccagctacgtggagcagatcgagaagaaacaaaaagtgtgaacatggctacgcctcttattcccgctccatgccgagccccgactgcagcattcaagccccacggaggtgagccgggc
acagaaacacagcagcgggagcagcaacaccagcactgccacagatctacacaaatgagctggaaaagaatcgagagctcatctgcgcctttgttgaacgcttaaaagtctgattccactaggaccagactgcacccggcacacaacttggtttgctcaacaagccaaagc
acacatcaagaacttgaagaagctgaagaaaaagccagcaccagctcgagaatttgaacgagaacagagatttttaaaagtggcgcactggaacagctgcagggctctcaggagatggaacgaatcgaatggacagcattggatcaactattctcagatcgttctgattcagagcga
gaggagattgaagtggatgttgaaagcagagatttcccatggagaagtggacaataatgaccaccagcatcagtgacattgatgaccacagcagcctgccgagatttggagtgacgaggggtactccagtgccagtgtaaaacttctcattcactcatag