



Reference: PTXBC035212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Flt4-L; LMPH1D; VRP; vascular endothelial growth factor C; FLT4 ligand DHM; vascular endothelial growth factor-related protein		

For research use only.

GENE DESCRIPTION	vascular endothelial growth factor C
GENE ID	7424
ORIGIN SPECIES	Human
ACCESSIONS	BC035212
VECTOR	pDONR223
NCBI SYMBOL	VEGFC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]