



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM171-transmembrane protein 171 Gene

Reference: PTXBC035310

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRP2; transmembrane protein 171; proline-rich protein PRP2		
For research use only.			

GENE DESCRIPTION	transmembrane protein 171
GENE ID	134285
ORIGIN SPECIES	Human
ACCESSIONS	BC035310
VECTOR	pDONR223
NCBI SYMBOL	TMEM171
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]