

KCNE1L-KCNE1-like Gene

Reference: PTXBC035330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KCNE1L; potassium voltage-gated channel subfamily E regulatory beta subunit 5; AMME syndrome candidate gene 2 protein; AMMECR2 protein; KCNE1-like; cardiac voltage-gated potassium channel accessory subunit 5; potassium channel subunit beta MiRP4; potassiu		

For research use only.

GENE DESCRIPTION	KCNE1-like
GENE ID	23630
ORIGIN SPECIES	Human
ACCESSIONS	BC035330
VECTOR	pDONR223
NCBI SYMBOL	KCNE1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaactgcagcgagagccagcggctgcgaacctctctgagccgctgtgtctcgcagctgcaccacggggaatgccagcggcttggcgctggccctctgccagcatgggcatgggggtcgtgaccttctgtggcgccgaggtgaccagcgccaagggcgacgacgc
ctatctctacatctgctcatcatgatctctacgcctgttggccggagggcctcatctgtgacctacacccgctcccgtaagctcgtcagggccaaggagagccgtccaggttgcgcgagcagcaatgggccccggaggcgccctgaccgccgacgccgaggtgcccggggc
tccagggcgagggcgccgagcgttgctccgaggggctgctgcctcgcacggcgctgagcgggtctaa