

TMEM139-transmembrane protein 139 Gene

Reference: PTXBC035517

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	transmembrane protein 139		
For research use only.			

GENE DESCRIPTION	transmembrane protein 139
GENE ID	135932
ORIGIN SPECIES	Human
ACCESSIONS	BC035517
VECTOR	pDONR223
NCBI SYMBOL	TMEM139
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgccaatgcacttactggggagactggagaagccgttctctctctgtgctgcgcctctctactggggctgcttctgctggcataaagacggacatcaccccggtgtattttcttcacattgggtggcttcttctgttgcctatctctggtccggtttctggaatggggcttc
gggtccagctccaatcaatgcagactgagagccagggccctcaggcaatgcacgggacaatgaagccttgaagtggcagctctatgaagagggcgtgggtgggactagaatccagtgccgcccccaagagttggaccaaccacccccctacagcactgtgtgataccccagcacc
tgaggaggaaacacctagccatccagaggggtccaggagagccaaactggaacagagcgcaatggcctcagaggggtccatggcccaggaaggagccctggaagagctccaatcaaccttcggcttcggggaccacgggctgtgtccactgctctgatctgcagagcttggcgg
cagtcgccacattagagcctctgactccacccctgcctatgatgtctgttggcaccctgatgatgatgtgtttttatgaggacaactggcaccctctaa