

TCEAL8-transcription elongation factor A (SII)-like 8 Gene

Reference: PTXBC035573

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WEX3; transcription elongation factor A protein-like 8; TCEA-like protein 8; transcription elongation factor A (SII)-like 8; transcription elongation factor S-II protein-like 8; transcription elongation factor A like 8		

For research use only.

GENE DESCRIPTION	transcription elongation factor A (SII)-like 8
GENE ID	90843
ORIGIN SPECIES	Human
ACCESSIONS	BC035573
VECTOR	pDONR223
NCBI SYMBOL	TCEAL8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaaagtcttgaagaaatgagggaaaaccacagaacatgccaaaggccgagggaagatgccctttggaggatgtaccacaggaggcagaaggaaatcctcaacttccgaagaaggcgttaagccagaagcagaaggaaacccagaggaggccgaatcagcctggc
cagggaattaaaggacacacccgttaggcatttggaccctgaagaaatgataaggagtagatgagcttgaaaggcttagggagagataagaagagtaagaacaagtgtgtgatgacattggaagcaagacattcacgcagccgtcttatcctgtgtgcttaggccttga