



Reference: PTXBC035588

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP2B; DAR; sarcoplasmic/endoplasmic reticulum calcium ATPase 2; ATPase Ca++ transporting cardiac muscle slow twitch 2; ATPase, Ca++ dependent, slow-twitch, cardiac muscle-2; SR Ca(2+)-ATPase 2; calcium pump 2; calcium-transporting ATPase sarcoplasmic ret		

GENE DESCRIPTION	ATPase, Ca ⁺⁺ transporting, cardiac muscle, slow twitch 2
GENE ID	488
ORIGIN SPECIES	Human
ACCESSIONS	BC035588
VECTOR	pDONR223
NCBI SYMBOL	ATP2A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1