

WWC3-WWC family member 3 Gene

Reference: PTXBC035601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein WWC3; BM042; WWC family member 3		
For research use only.			

GENE DESCRIPTION	WWC family member 3
GENE ID	55841
ORIGIN SPECIES	Human
ACCESSIONS	BC035601
VECTOR	pDONR223
NCBI SYMBOL	WWC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagtgccgtgcttccgtatcagtcgaaggagccctcgtgtttgccacctttaccttgaacctccccctgcctccctgctgtgtccgcttttgcagctcaatgcagccatgacaaggaaagaaaaagacaagggaagccagagagccgcgcagttctctgcaggtgcagatgcaggcag
tggagggtggcctgagcaggcagaaggacaccaagcgcctatgttgccttgcattcatgacgtgtgttggagcttctgactagttcagactgccacgcccaacccagaaaaatacccatgccagaaaagtgaagtcctaggtgtttccatctatgtttcaatctgtccatctaccaggcctc
gcgataaaaaacaaaaaacgcctgccaggttttagaagcagttctgtgtctcaaacatcaggatcctgccaccaggggtctttgaaatagtagcacatgtaaaagggaatttggctttcacttcataatcactga