



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
 contact@proteogenix.fr
 www.ProteoGenix.science

FFAR3-free fatty acid receptor 3 Gene

Reference: PTXBC035657

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FFA3R; GPR41; GPR42; free fatty acid receptor 3; G-protein coupled receptor 41		
For research use only.			

GENE DESCRIPTION	free fatty acid receptor 3
GENE ID	2865
ORIGIN SPECIES	Human
ACCESSIONS	BC035657
VECTOR	pDONR223
NCBI SYMBOL	FFAR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr