

FAM72A-family with sequence similarity 72, member A Gene

Reference: PTXBC035696

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM72A; LMPIP; Ugene; LMP1-induced protein; latent membrane protein 1-induced protein; family with sequence similarity 72 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 72, member A
GENE ID	729533
ORIGIN SPECIES	Human
ACCESSIONS	BC035696
VECTOR	pDONR223
NCBI SYMBOL	FAM72A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctaccacattttagttcaaggacaggtgctgtccatcctgtgttgcaaatctgttaacaagtgtcagcttaggggaatgaaggctgttttctggctgatactgaaatagacctttctctacagacatccctctaccaacgcagtgacttcaactggaagatgctatccacaa
atctgcaaatgtaactgaaggacatgcatgttttaaatgtgggaacattgttaggtatcatgtgattgttcctgtagttcctgtcttcttctgcaacaacggacacttctggatgtttcacagccagcagtttatgatattaacagactagactccacaggtgtaaacgtcctactttgggca
actgccagagatagaagaggtacagatgaagatgtgttaaatatctcagcagaggagtgattagataa