

C2orf80-chromosome 2 open reading frame 80 Gene

Reference: PTXBC035737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C2orf80; GONDA1; gonad development associated 1; chromosome 2 open reading frame 80		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 80
GENE ID	389073
ORIGIN SPECIES	Human
ACCESSIONS	BC035737
VECTOR	pDONR223
NCBI SYMBOL	C2orf80
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagaaggctcataaagaaggaatgaaaaagctctgggagattatattggcatcagactcgggaaatgaattgacccaaaaggaagacggcaactcactttctagatgatatggcacacatgacttggccatcagtggttcttcaatggctggatccctcagaagacttaa
cttggctggagtgagggaactgaaaataccactccatggcagaccatataatcgaatcgtagagaacgagaagctatgattttatcatcttctgctggaatcttaataacagatcccgattgagggaagctcttaaaattatgggctgattcttctgcccgaattctgtaccatcaagggtccc
cgagtttcatctctccgctctcttgcaccccttggccatgtaacagacccaaagcagcagcatagcccgcaaacagagtgtaagtaagaaggttaacaacaataaaaaatgccaccagcatctctgcaagggaagcaaatgccacagaatggaaatcatcacaagggttttcaga
cacacagccaagcacaagtcacttag