

LOC149643-hypothetical LOC149643 Gene

Reference: PTXBC035742

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf227; spermatogenesis-associated protein 45; UPF0732 protein C1orf227; spermatogenesis associated 45		

For research use only.

GENE DESCRIPTION	hypothetical LOC149643
GENE ID	149643
ORIGIN SPECIES	Human
ACCESSIONS	BC035742
VECTOR	pDONR223
NCBI SYMBOL	LOC149643
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatctataaacagaacattgaaataatgaaaaacatggagtaagcaaacatctctggaggagataaacaacaaagcgtgaatccaactgcttggtggaacgaagcaatcaagtcagcttactgagagttcaaaagaggcacttccggatgcctatcagtcctttactgatacc
acaaccaaagagcctgttcccaacagtggcaggagctcttgatcaagctgagtcctctgctcacatggagagaaagcactttccacaaaaataatgccatatttggataa