

TM4SF20-transmembrane 4 L six family member 20 Gene

Reference: PTXBC035754

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO994; SLI5; TCCE518; transmembrane 4 L6 family member 20; transmembrane 4 L six family member 20		

For research use only.

GENE DESCRIPTION	transmembrane 4 L six family member 20
GENE ID	79853
ORIGIN SPECIES	Human
ACCESSIONS	BC035754
VECTOR	pDONR223
NCBI SYMBOL	TM4SF20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgctgcaagatggacatcctgcaatggattcagctgctgttctactgctgttaggagtagttctcaatgtgatacctctaattgtcagcttagttgaggaagaccaattttctaaaacccatcttctgtttagtggtggtgccaggaattataggagcaggtctgatggccatt
ccagcaacaacaatgtccttgacagcaagaaaaagagcgtgctgcaacaacagaactgggaatgtttcttcatcattttcagtgatcacagtcattggtgctctgtattgcatgctgataccatccaggctctcttaaaaggctcctcatgtgtaattctccaagcaacagtaattgccaattgt
gaattttcattgaaaaacatcagtgacattcatccagaatcctcaactgacgtggttttcaatgactctgtgacactcctactgggttcaataaaccaccagtaacgacaccatggcgagtggtggagagcatctagtctcacttcgattctgagaagaaacacataggcttatccactct
cagtatttttaggtctattgctgttggaattctggaggtcctgtttgggctcagtcagatagtcacgtttccttggctgtctgtgtggagtccttaagcgaagaagtcaaatgtgtag