



Reference: PTXBC035784

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TM14; fibulin-7; FIBL-7; fibulin 7		
For research use only.			

GENE DESCRIPTION	fibulin 7
GENE ID	129804
ORIGIN SPECIES	Human
ACCESSIONS	BC035784
VECTOR	pDONR223
NCBI SYMBOL	FBLN7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]