

C19orf66-chromosome 19 open reading frame 66 Gene

Reference: PTXBC035817

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0515 protein C19orf66; RyDEN; repressor of yield of DENV protein; Repressor of yield of Dengue virus; repressor of yield of DENV; chromosome 19 open reading frame 66		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 66
GENE ID	55337
ORIGIN SPECIES	Human
ACCESSIONS	BC035817
VECTOR	pDONR223
NCBI SYMBOL	C19orf66
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcaggaaggtgtggagctggagaagagcgtccggcgccctcgggagagaagtttcatgggaaggtatcctccaagaagggcggggcctctgatgaggaattcggcagcgaccacacgggagtggggcgtccatcgtgtacggggtaaagcaaaaagatggccaagaactaag
taacgatctggatgccaggaatccaccagaagatatgaagcaggaccgggacattcaggcagtgggcagctcctcctcctgacagagaagccaacctacgcattttcaacgtgccaggagacattatccctgctgtggaccggcagtttgcctgctcctctgcgaccacgtctgg
tggcgccgctgcccagcggaaggaggtatccgggtgcccgaatgccggaagcgtacgagccagtgccagctgacaaagatgtgggctgctgagttccactgccgaagtgtcggcacaactccgcacccacactcactcctgctcagctgcccactgctacaaccggcga
gagccccacgtgcctgggacatcctgtgtctaccccaagagccggaagcagaaccactgcccagaagtgtcctcccccagcaacctacattagcagtggtccactgtggccactgcttgagccagggtggcctcctggaagacctggacaacctatcctggaggactgaagg
aggaggaggaggaagaggagggtggaggacgaggaggcggggccaggagtgga