

ERAF-erythroid associated factor Gene

Reference: PTXBC035842

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERAF; EDRF; alpha-hemoglobin-stabilizing protein; alpha hemoglobin stabilising protein; erythroid differentiation associated factor; erythroid differentiation-related factor; alpha hemoglobin stabilizing protein		

For research use only.

GENE DESCRIPTION	erythroid associated factor
GENE ID	51327
ORIGIN SPECIES	Human
ACCESSIONS	BC035842
VECTOR	pDONR223
NCBI SYMBOL	ERAF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcttctaaggccaataaggatctcattccgaggtgaaggagttcagcgttctgtgaatcagcaggtcttcaatgacctctcgtctctgaagaagacatggtgactgtggtggaggactggatgaacttctacatcaactattacaggcagcaggtgacagggagccccaag
agcgagacaaggctctgcagagcttcggcaagagctgaacactctggccaaccttctctggccaagtacagggacttctgaagtctcatgagctcccgagtcacccaccgcccctctctag