



Reference: PTXBC035893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FPN1; HFE4; IREG1; MST079; MSTP079; MTP1; SLC11A3; solute carrier family 40 member 1; iron regulated gene 1; solute carrier family 11 (proton-coupled divalent metal ion transporters), member 3; solute carrier family 40 (iron-regulated transporter), member		

For research use only.

GENE DESCRIPTION	solute carrier family 40 (iron-regulated transporter), member 1
GENE ID	30061
ORIGIN SPECIES	Human
ACCESSIONS	BC035893
VECTOR	pDONR223
NCBI SYMBOL	SLC40A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]