



Reference: PTXBC035900

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sodium-independent sulfate anion transporter; solute carrier family 26 (anion exchanger), member 11; solute carrier family 26 member 11		

For research use only.

GENE DESCRIPTION	solute carrier family 26, member 11
GENE ID	284129
ORIGIN SPECIES	Human
ACCESSIONS	BC035900
VECTOR	pDONR223
NCBI SYMBOL	SLC26A11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]