

DDX55-DEAD (Asp-Glu-Ala-Asp) box polypeptide 55 Gene

Reference: PTXBC035911

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent RNA helicase DDX55; DEAD (Asp-Glu-Ala-Asp) box polypeptide 55; DEAD box protein 55; DEAD-box helicase 55		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 55
GENE ID	57696
ORIGIN SPECIES	Human
ACCESSIONS	BC035911
VECTOR	pDONR223
NCBI SYMBOL	DDX55
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagccccagagaacacacagcggacctctgccaaaactcaagtccatggccctggctgacagagctgtgttgaagggcatgaaagctttgtcatatgtccaagcttatgcaagcatgaatgaacctgatttcagattaagatcttgattgccagccttctcgagtttt
gcccctctgaggatccccagatccagaattgagagggaagcagttccagattttgtcccgtggacgttaataccgacacgattccattaaagataaaatcagagaaaagcagagcagaaaactctggagcaacaaaagagagaaaacagaaaatgaaggagagaagaaatt
cataaaaaataaagcttggtaaacgagaagcacaagaaagaaaagaaaatgaatgagaaaagaaaagggagaggggtctgatattgaagatgaggacatggaagaactcttaatgacacaagaactcttgaaaaaactaaagaaagcacaataactgaagaagaatttg
agaaggcctgtgacaactggcaaaagaacaatcaagacagtgatttagggatctcagattggaagatgactgctaa