

HCST-hematopoietic cell signal transducer Gene

Reference: PTXBC035931

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP10; PIK3AP; hematopoietic cell signal transducer; DNAX-activation protein 10; kinase assoc pro of ~10kDa; kinase assoc protein; membrane protein DAP10; phosphoinositide-3-kinase adaptor protein; transmembrane adapter protein KAP10		

For research use only.

GENE DESCRIPTION	hematopoietic cell signal transducer
GENE ID	10870
ORIGIN SPECIES	Human
ACCESSIONS	BC035931
VECTOR	pDONR223
NCBI SYMBOL	HCST
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccatctgggtcacatcctctctctgtttgtccagtggtgcagctcagacgactccaggagagagatcatcctcctgcttttacctggcacttcaggctctgttccggatgtgggtccctctctgcccgtcctggcaggcctcgtggctgctgatgcggtggcatcgtgc
tcacgtggggcggtgttctctgtgcagcggccacgcgcagccccgcccaagaagatggcaagtctacatcaacatgccaggcaggggctga