



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SERF1A-small EDRK-rich factor 1A (telomeric) Gene

Reference: PTXBC035932

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4F5; FAM2A; H4F5; SERF1; SMAM1; small EDRK-rich factor 1; SMA modifier 1; protein 4F5; small EDRK-rich factor 1A (telomeric); spinal muscular atrophy-related gene H4F5; small EDRK-rich factor 1A		

For research use only.

GENE DESCRIPTION	small EDRK-rich factor 1A (telomeric)
GENE ID	8293
ORIGIN SPECIES	Human
ACCESSIONS	BC035932
VECTOR	pDONR223
NCBI SYMBOL	SERF1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgctgggaatacaacgagaacttgcccgccagaaaaatgaagaaaccaggaaattgcaagggaagaggaaaggagtagcttgactgccttcagagaaagcagagtctggggccgaaaactgagagcaagatgtcagctggggccacacctccctctgaaggct
ccaagggaagaactctgtcttcctccagctgctgtggtccagctatcttgcttatgccagacataactcctatctcctccttctcttctgtgtctctttctctctctccctcttttatagagacatttctctgttgatttag