



Reference: PTXBC035960

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLEPP1; NPHS6; PTP-OC; PTP-U2; PTPROT; PTPU2; R-PTP-O; receptor-type tyrosine-protein phosphatase O; PTP phi; PTPase U2; glomerular epithelial protein 1; osteoclastic transmembrane protein-tyrosine phosphatase; phosphotyrosine phosphatase U2; protein tyro		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, receptor type, O
GENE ID	5800
ORIGIN SPECIES	Human
ACCESSIONS	BC035960
VECTOR	pDONR223
NCBI SYMBOL	PTPRO
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr