

## UMOD-uromodulin Gene

Reference: PTXBC035975

|                |                                                                                                         |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ADMCKD2; FJHN; HNFJ; HNFJ1; MCKD2; THGP; THP; uromodulin; Tamm-Horsfall urinary glycoprotein; uromucoid |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | uromodulin                                                                                       |
| <b>GENE ID</b>          | 7369                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC035975                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | UMOD                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggggcagccatctctgacttgatgctgatgggtgggtggcctctgtgttcacacactgcagccactgacacctcagaagcaagatggtgctctgaatgtcacagcaatgccacctgcacggaggatgaggccgttacgactgacacctgtcaggagggtcacggcgatggc  
ctgacctgctgggacctggatgagtgccgcatctctggagctcacaactgctccgcaacagcagctgctaaacacgcccaggtctctctctgctgctgccccgaaggcttccgctgtgccccggtctcggctgcacagacgtggatgagtgctgagcctgggcttagccactgc  
cacgacctggccacatgtgtcaatgtgtgggtgggagctactgtgctatgccccggtgctacgggggagtgatggcactgtgagtgctctccgggctctgccccgggggtggactgctgccccgaggggcagcgctcgtgtgctggatccgtgtcaggcgaccgcacc  
ctggacgagtactggcgagcaccgagtgacggggagggtctacgctgctgcacagcggacctgcgctgggtgacgcccccatcgtccagcgacgaggggcatgctgagccgcaaggcctgcgcgactggagcgccgactgctgctgtgggatgctcctccagggtgaaggcc  
tggccggcgctactacgtctacaacctgacagcgccccgagtgacctggcgtactgcacagccccagctccgtggaggggagctgtgaggagtgagtagacgaggactgcaaatggaataatggcagatggcactgccaatgcaaacaggactcaacatcactgatat  
ctccctctggagcacagcgtggaatgtgggccaatgacatgaaggtgtcgtctggcgaagtgcacgtgaagagtgctggcttcgacaaggtctctcatgtacctgagtgacagccggtgctcgggcttcaatgacagagacaaccgggactgggtgtctgtatgacccccagccggg  
atggccccctgtgggacaggtgtgacgagggaatgaaacccatgccacttacagcaacacctctactctggcagatgagatcatcctcgtgacctcaacataaatacaacttggcatgctctacccccctggacatgaaatgacgctgagaccgccccacagccaatgtcagtgctcta  
aacatcagatggggcgccggcgtgttaccgtgcggatggcgtctctcagacccttctacacgacccctaccaaggctcctcgtgacactgtccactgaggtttttctctacgtgggaccatgttggatggggggcagactgtcccgatttgactgctcatgaccaactgctat  
ggcacaccagtagcaatgccacggacccccgaagtacttcatcatccaggacagatgccacacactagagactcaactatcaagtggtggagaatggggagtgctccagggccgattttccgtccagatgttccggttggctggaaactatgacctagtctacctgacgtgaaat  
ctatctctgtgacacatgaatgaaaagtgcagccctacctgctctgggaccagattccgaagtgggagtgctatagatcaatccggtgctctgaaactgggtccatcacacggaagggtgtccaggccacagtgctcaaggccttttagcagcttggggctctgaaagtctgctgctctg  
cttctctggccacctgacctgacttttcagtga