

HSD17B11-hydroxysteroid (17-beta) dehydrogenase 11 Gene

Reference: PTXBC036001

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	17-BETA-HSD11; 17-BETA-HSDXI; 17BHSD11; DHRS8; PAN1B; RETSDR2; SDR16C2; estradiol 17-beta-dehydrogenase 11; 17-beta-hydroxysteroid dehydrogenase type XI; CTCL tumor antigen HD-CL-03; CTCL-associated antigen HD-CL-03; cutaneous T-cell lymphoma-associated a		

For research use only.

GENE DESCRIPTION	hydroxysteroid (17-beta) dehydrogenase 11
GENE ID	51170
ORIGIN SPECIES	Human
ACCESSIONS	BC036001
VECTOR	pDONR223
NCBI SYMBOL	HSD17B11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaattctctggacatcctctgcttccctgtaactgacgtctgtctccctagagtcctctgtgaagctttttattcctaagaggagaaatcagtcaccggcgaaatcgtgctgattacaggagctgggcatggaattgggagactgactgcctatgaattgctaaactaaaagcaagct
gggtctctgggatataaataagcatggactggaggaacagctgccaaatgcaaggagctgggtgccaagggttcacacctttgtgtgtagactgcagcaaccgagaagatattacagctcgcgaagaagggtgaaggcagaattggagatgttagtatttagtaataatgctggtgtagtc
tatacatcagattgtttgctacacaagatcctcagattgaaaagacttttgaggttaagtactgacatttctggactacaaaggcatttctctgcaatgacgaagaataaccatggccatattgtcactgtggtctcggcagctggacatgttccgtccctcttactggttactgttcaagc
aagtttgctgctgttgattcataaaactttgacagatgaactggctgacctacaataactggagtcaaaacaacatgtctgtcttaattcgtaaacactggcttcacaaaatccaagtaacaaatttggaccactctggaacctgagggaagtgttaaacaggctgatgcattgggttct
gactgagcagaagatgattttattccatctcttatagctttttaacaacattggaaggatcctctgagcgtttctgagcttttaaaacgaaaaatcagtttaagtttgatgcagttattggataaaaatgaagcgcaataa