

## TEX11-testis expressed 11 Gene

Reference: PTXBC036016

|                |                                                                                                              |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SPGFX2; TGC1; TSGA3; testis-expressed sequence 11 protein; testis expressed sequence 11; testis expressed 11 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | testis expressed 11                                                                              |
| <b>GENE ID</b>          | 56159                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC036016                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TEX11                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggactttaagaagttgttgaacctggttacaatgataatcacataccagaggaagcgaatgatgactctcagcgacatagcaaatcaacagggagctctatggctgaataacagacattcagattgaagaatggcagtaaacctatggaactgggcacttaccataggag  
gaggttggctgttaaatgaagagcagaaaattgattacattatgttctgcaagttgctgagttgttgaagcctcttggctcagaacaaagtattcaacgactgattatgatgaatatgagaatggaagaatggttggatgctggaattttctaatcctgatgaatgtttcaagctg  
ctgtggcagctctggagcaattatcgtcaaatattcaagagctccctgaggtgactgacattgagagattactgttgaagttgaccacttcagagtgcttcttaccagcagagtcagcagttgctcaaggggattttcaagagcattctgtgtgactgcaatgttaagat  
gttgaaggctccccagatgacttcaagcttcatctctgtttacaacttttgagtagaaccagaagaataataatgaagaaagtcttctgcttagcgaagctatgatttgggaatggataagaatctactgggcccagaatgctggctaaagtctacggctattagcc  
acgaattatttggatgggatgacaccaatattatgataaggctctcaatgctgttaaacctagcaaacagggaacatttaagttctctgggctttcttaaaaatgaaaatcctcttgaaggcgaacatctaatgaagaactccttgaagctgtatggaatactactctgacatgcctta  
gacttctgtcgaacattgctaaactgctgatgcatgaaagagaatctgttgggtttcatttctgacgattatcatgaacgttttaagtcacggaataattggaagaattctgatactccatactgacatgcttttacaaggaaggaagaactcttgcgaaggagaattgaagaatct  
tttttagctcaccacaacggaagacaactgacagcagaatcaatgaactggttacacaacattctgtggagacaagctgccagtagtttgggtgacaaaattacactgatgcctacaatggtgactatttctctgaggttttattcaactgataaaatggaatctggacttcaccaagctgcagag  
gaacatggctgtcttactgaatttgcacaactgataaggccaaagaggcagtgggcagaagctgaacgacatgacctaggaacgttttactcaatttatataatcaagattgacgtcatagagggcaactctgaagagcttgcaggcgaataattcttagagataatataacagat  
gaagagtcagaagataatgatctagttgcagagagaggttcacctaccatgcttcaagtttagctgccagtttgcctctagagaatggacaacaaattgtgacagaaaagcttggaaatatttagctcaacattcagaagaccagggaacaggttctacagctgtaaggtgttctgcttctt  
cttccaaaattgctgaaatgccggaatctgaagataagaagaagaagaatggatcgaattttagcttgcctgaatagagccttggtaaaccttctcagccttgggtgaagaagccttaagtttggagtcagagcctaatgaagctcagtggtttcgaagaaacagcttggacttggctgtgcaat  
gtgacaagatcagtgatgagagagtttttatacttttataagatgtcccaagtttgcctctctgacgaagtaattctgattgcacgaaaacatgtttactatgacagttgcagttgacgttagagcaaggggaaaaagcttcaacagcctttgaacagaccatgttctgagctgcact  
gaggagatccagacatgcaatgacatccataatttctgaaacaaacaggggaccttctcaaatgattcatgtgagaattgcttctgctgacgagtttgaagttgagccaaatgaaatgatccattactggaagcttctggaatcagtggtggagttgctcattagaactaaacatttga  
aacaattgcaataatagcaatggaagcctgcacactatcctttagttgctctcaaggccttgaagaaaggcttatttctctacaaaagggaagaaacaaatgatatatcacaaatgacagcaaatgtatgcacaacttggtaacctctcagtgccagatggggcgtcgaatgttagagctgtcc  
cctggaagaagtttgggctattttgaagacgtctgagccacattagccgactaaagactaccgaagaatggagattctctggctgagtgacagtcctggaatccggagtaactattttagcaggagcaagatgcatctgctgaaaggtgtgtggcctggccttgccttcttaacc  
accttacctcttcaaggaaagctatgaaactcagatgaatgctgtatagctcagctgtgtgaagcattgagtaacaacaaaggcccagttttcatgaacatggctactggagcaagtcagattag

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)