



Reference: PTXBC036037

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	alpha-(1,3)-fucosyltransferase 11; alpha (1,3) fucosyltransferase; fuc-TXI; fucT-XI; fucosyltransferase XI; galactoside 3-L-fucosyltransferase 11; fucosyltransferase 11		

For research use only.

GENE DESCRIPTION	fucosyltransferase 11 (alpha (1,3) fucosyltransferase)
GENE ID	170384
ORIGIN SPECIES	Human
ACCESSIONS	BC036037
VECTOR	pDONR223
NCBI SYMBOL	FUT11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]