

EIF5A2-eukaryotic translation initiation factor 5A2 Gene

Reference: PTXBC036072

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIF-5A2; eIF5AII; eukaryotic translation initiation factor 5A-2; eIF-5A-2; eukaryotic initiation factor 5A; eukaryotic translation initiation factor 5A2		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 5A2
GENE ID	56648
ORIGIN SPECIES	Human
ACCESSIONS	BC036072
VECTOR	pDONR223
NCBI SYMBOL	EIF5A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagacgaaattgattcactactggagatgccggggctccagcactaccctatgcagtgctcgcccttgcgcaaaacggcttcgtggtgctcaaggacgaccatgcaaaatagtgagatgtcaactccaaactggaagcatggtcatgccaaggtcacctgttgaatt
gatatttcacgggcaaaaaatagagatattgtccttctactcacaacatggatgttccaaatattaagagaaatgattatcaactgatgcattcaagatggttaccttccctgctgacagaaactggtgaagttcgtgaggatcttaactgccagaaggtgaactaggcaagaaataga
gggaaatacaatgcaggtgaagatgtacaggtgtctgtcatgtgtgcaatgagtgagaatatgctgtagccataaaacctgcaataa