

GAGE1-G antigen 1 Gene

Reference: PTXBC036094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT4.1; GAGE-1; G antigen 1; MZ2-F antigen; cancer/testis antigen 4.1; cancer/testis antigen family 4, member 1		

For research use only.

GENE DESCRIPTION	G antigen 1
GENE ID	2543
ORIGIN SPECIES	Human
ACCESSIONS	BC036094
VECTOR	pDONR223
NCBI SYMBOL	GAGE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgaggagatcgacctattattggcctagaccaaggcgctgtacagcctctgaaatgattggcctatgcggccgagcagttcagtgatgaagtggaccagcaacacctgaagaaggggaaccagcaactcaacgtcaggatctgcagctgctcaggaggag
aggatgaggagcactctgcaggtcaaggccgaagcctgaagctgtagccagggaacagggtcaccacagactgggtgtgagtgtagatggtcctgatgggcaggagatggaccgccaaatccagaggaggtgaaaacgcctgaagaagatgagggtcactatgttgc
cagactgggattctctggcttttaatgaacaattgcttcttaaatcttccccacggaaccttga