



Reference: PTXBC036187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	160-KD; POP101; SRM160; serine/arginine repetitive matrix protein 1; SR-related nuclear matrix protein of 160 kDa; Ser/Arg-related nuclear matrix protein; serine/arginine repetitive matrix 1; serine and arginine repetitive matrix 1		

For research use only.

GENE DESCRIPTION	serine/arginine repetitive matrix 1
GENE ID	10250
ORIGIN SPECIES	Human
ACCESSIONS	BC036187
VECTOR	pDONR223
NCBI SYMBOL	SRRM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr