

DLX1-distal-less homeobox 1 Gene

Reference: PTXBC036189

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein DLX-1; distal-less homeo box 1; distal-less homeobox 1		
For research use only.			

GENE DESCRIPTION	distal-less homeobox 1
GENE ID	1745
ORIGIN SPECIES	Human
ACCESSIONS	BC036189
VECTOR	pDONR223
NCBI SYMBOL	DLX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttggcgcccaaccagcaaatgtctctctccatgtccacgggcactactcactgactgtttacactggcgccatcgcagcccgacggcgccctacagctcagcctctctctccgaccgctgggtaccctacgtcaactcggtcagcagccacgcatcag
ccctacatcagttcgggtcagtcctaccggcgagcgcagcctcgcagagccgctggaggaacccaggggcgactcggagaagagcagcgtgtggaaggcggtgaagtgcgcttcaatggcaagggaagaaagatccgtaaacccaggacgattattgcagttgcagtt
gcagccttgaaccggaggtccagcaaacctcagctacgtctgtccggagagggcgagctcgcggcctctttgggactcacacagactcagggtcaagatctgttccaaaacagcgatccaagttcaagaagctgatgaagcaggtggggcggtcctggagggtagtcggtt
gccaacggcggcgccctgtctgtgctcccccacccgtgcccggcgctggaaccctaactcttcatccgggaagggtcaggaggaacgcgggctctatatcccgactacacatcgtgtacccttcagcgaccaagaagctatgcagcaacccaactatgtga