



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MGC39545-hypothetical protein LOC403312 Gene

Reference: PTXBC036197

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

putative uncharacterized protein MGC39545

For research use only.

GENE DESCRIPTION	hypothetical protein LOC403312
GENE ID	403312
ORIGIN SPECIES	Human
ACCESSIONS	BC036197
VECTOR	pDONR223
NCBI SYMBOL	MGC39545
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1