

PLAC1L-placenta-specific 1-like Gene

Reference: PTXBC036256

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLAC1L; TMEM122; oocyte-secreted protein 2; placenta-specific 1-like protein; testicular tissue protein Li 142; transmembrane protein 122; oocyte secreted protein 2		

For research use only.

GENE DESCRIPTION	placenta-specific 1-like
GENE ID	219990
ORIGIN SPECIES	Human
ACCESSIONS	BC036256
VECTOR	pDONR223
NCBI SYMBOL	PLAC1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtagaagcttctgctcctcgtgtcttgattggaccgggtgctgagaaacctccatgtgaaataagttgctctctgactggttgatggtctcagttatccagttgcagaaagcagaaatctgtatatatttgcggatgaattacatctgggaatgggctgccctgcaaatcggataca
tacatatgtatagatttatatactgttcgtgattgtggcatcaggacaagggtagttctgaggaaactctcttttcaaacaggagctgtactttacccaaggaatagatcatgacctcaggaaatccatttggaggtttccacctctaggaaatcagttgtggttacaccagtttctactga
gaatgaaataaaattggatcctagtcctttattgctgacttcagacaacagcagaagagttaggattattatctctagtcctcaaaactgctctga