

TMEM192-transmembrane protein 192 Gene

Reference: PTXBC036301

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 192		
For research use only.			

GENE DESCRIPTION	transmembrane protein 192
GENE ID	201931
ORIGIN SPECIES	Human
ACCESSIONS	BC036301
VECTOR	pDONR223
NCBI SYMBOL	TMEM192
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgggggcaggatggaggacgttccttgatcacccagagtattgaagacgacccactctggatgccagcttctccacaccactattacaagctcactttagaccccgattccatctctctacagtcacatagtgaaattctgtgtttattcatctcgtgtttgttt
ttagcatttttaacagggtgtgtgtttcttctctaatccaaatgaggacaagtgccaggaaattacacaacccattgaaagttcagacggttataatcttgggaaagtattttgtggattctccatttactcttgaatgctacatccagatcacccacagaaaatcagaaccgaggctata
acttgatctaccgatcaacaaggcatctcaagagacttgcgttgatgatacagctcttggcaacacagtgcttctctcactactgtgcatgcagcactctctccagagcctggcagattgtatcttgacctattctggccatcttggcactggactcaatctgttccctgataatgtctcctcattta
cacagtgaataacgggagatttaataaagcctaaccagagcctgatacttgaagaagaaaaatctatgcttaccagcaatattacctggagactggattcagaactattcaagcctagaagaattgttgaagaagcaggagacaccattgaatacctgaagcgacacaatgcgct
gctgagtaagcgattgttgctctcacttctcagacctgggctgtcagccaagtagaacgtga