



Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TLR4 interactor with leucine rich repeats		
For research use only.			

GENE DESCRIPTION	KIAA0644 gene product
GENE ID	9865
ORIGIN SPECIES	Human
ACCESSIONS	BC036337
VECTOR	pDONR223
NCBI SYMBOL	KIAA0644
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]