

C10orf81-chromosome 10 open reading frame 81 Gene

Reference: PTXBC036365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PH domain-containing protein C10orf81; C10orf81; HEL185; bA211N11.2; pleckstrin homology domain-containing family S member 1; PH domain-containing family S member 1; epididymis luminal protein 185; pleckstrin homology domain containing S1		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 81
GENE ID	79949
ORIGIN SPECIES	Human
ACCESSIONS	BC036365
VECTOR	pDONR223
NCBI SYMBOL	C10orf81
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacaaagtctccaggttaggcaaacacacacagattatcagaagccactcaagatgtgaaggagagaatcattcttactctcgaagtgttcttttagagttggataatattcgtccagtgattctgtggaatccattgaaactgatgtccagaccaggtctctggaag
aattgagtgctattatgagccaatggaatcctatttttcaagagacatccatgagctgtgtgagtagcagcaaaaggaggaacccagaccctccagagaccaggtggggacctccactgcaagaacaaggctcagggaattgattgtgtcttccctgccgatgtggaagcacagac
cacaaatgacaaaagggtgctgacctactaactgtgtgcaattgtctatattatcaataatatacccgatgaaagccaagtggagaaactgaacgttttcttctctctgtatcatcaactatctgtctctacagaagccacaggacgatatgtgtctcagtggaagggcccg
acgtttggatgcataattttgccacggagatccttctgcaagtgaatgacctgaaacccagagcctggaggaggtctcctgtttcttaccgggtccatccagaaggagaaattaaagcttaccatggcaggatcccaattcagagacattccatgctgcacctgtatgtccctcaa
atgccaaagtgtgcaccttctcagctggataagcctagactgaacagagctcccaagagagtgccggccattaaaaagagccagcagaaggagccagggagtaa