



Reference: PTXBC036439

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-203; HPF12; HTF12; TF12; zinc finger protein 92; epididymis luminal protein 203; zinc finger protein HTF12		

For research use only.

GENE DESCRIPTION	zinc finger protein 92
GENE ID	168374
ORIGIN SPECIES	Human
ACCESSIONS	BC036439
VECTOR	pDONR223
NCBI SYMBOL	ZNF92
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]