

C5orf35-chromosome 5 open reading frame 35 Gene

Reference: PTXBC036528

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf35; SET domain-containing protein 9; SET domain containing 9		
For research use only.			

GENE DESCRIPTION	chromosome 5 open reading frame 35
GENE ID	133383
ORIGIN SPECIES	Human
ACCESSIONS	BC036528
VECTOR	pDONR223
NCBI SYMBOL	C5orf35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctggcgcgtctgctgcggggcctgtggcagcgtggcgcgttacaagtacgctctgtccctggatgcactgaacctaaagccacaacccgaggaccccgatatttccagaggaaatccaaagacaaagtattcagatgaagatgtcctaggaacattactgaaagtgttcca
ggctctattcttaaatgattcaataaacaatcagaatcttgactatgcttcagaatctgttaatacaaaatcaagacctactggcagttgaacatcaaggggtgaaactgcgtgaaaacagacatcaacagcaaaagtacctttaaccagaagaattcttacaagacttgggttcagtgt
tgcccaagcaactagctcattgatttctgctggaagggtgtcttctactaaaggatgtgtacccaaaggcgcagtcgtatctatgtatcctggtacagtatacagaagtatgagccgatcttttccagtcattggaaatccgtttatttttagatgcctggatgggtactcattgatgggaat
gacaaaggatatacaaaagtgtgtacagatctgcaatgggaggatcgactggccctttaaaatgagtgatagtagatggcctaaagtcagaaatcataccctctggcgtgtgggacagtatgtcaacaattgttccaatgacagagcagctaatgtctgttatcaggaattgatgtgcct
gcagtttccctatagaactgaagcagtagcttccaacattgctacagctatgacaacaaggccacttcgatgtgtgtctgtgcacttagggacatcaatcaaggagaagagctttttcaactactacacaattgtcagctaa