

PTPN20B-protein tyrosine phosphatase, non-receptor type 20B Gene

Reference: PTXBC036539

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTPN20B; CT126; PTPN20A; bA142I17.1; bA42B19.1; tyrosine-protein phosphatase non-receptor type 20; protein tyrosine phosphatase, non-receptor type 20A; protein tyrosine phosphatase, non-receptor type 20B; protein tyrosine phosphatase, non-receptor type 20		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 20B
GENE ID	26095
ORIGIN SPECIES	Human
ACCESSIONS	BC036539
VECTOR	pDONR223
NCBI SYMBOL	PTPN20B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcacctaggacatttagacagagcctgtaaacgattatgagggaatgactctgaagcagaagacttgaattcaggagactttgcctcatcaagtcaggaaaacacacctagatcaaagggttttgaataaaagttaattcagagaaggtaaaactttctctcggaattccac
ataatgattatgaggatgttttgaagagccttcagaaagtggcagtgatccagcatgtggacagccagaggcccttcagaagagacaggtggagcagtgaggatgaggaggctgcaggccatcacaggctctctccctctactttctgtacgcgcgcaaaattgttctgaaggagaa
ctagatcagttggctcagattcggccattaatcaattttcatgagcagacagccatcaaggattgtttgaaatccttgaagaaaaaacagcagcgtatga