

PLAC2-placenta-specific 2 (non-protein coding) Gene

Reference: PTXBC036545

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLAC2; LINC00036; NCRNA00036; tissue differentiation-inducing non-protein coding RNA		

For research use only.

GENE DESCRIPTION	placenta-specific 2 (non-protein coding)
GENE ID	257000
ORIGIN SPECIES	Human
ACCESSIONS	BC036545
VECTOR	pDONR223
NCBI SYMBOL	PLAC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacctggcctgaccagctgatgcacactgcttcagacactcctgctggagcccgatccctgacaaggacctaggacattttgctcctgccagcctatgggagggagccttgagccttcagctctgctgtgactttgaggtttgtctccctctggggccctgggtgcctg
tcttcagtggaagcactgtgccaccttggaagctccatgggagccagagggcatcgcaagaagagaagcacagaagggcgaggagagacactcagaggcacttccgctcttgccaggacatttccagccacatttggccaagccgtgccccctgctggagcactttca
accttctctgcagctccaatacacctgggattgcagtcctccaggaagtctctcagattccctcctccagccagagagcacctag