

ASB17-ankyrin repeat and SOCS box-containing 17 Gene

Reference: PTXBC036573

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Asb-17; ankyrin repeat and SOCS box protein 17; ankyrin repeat and SOCS box containing 17		

For research use only.

GENE DESCRIPTION	ankyrin repeat and SOCS box-containing 17
GENE ID	127247
ORIGIN SPECIES	Human
ACCESSIONS	BC036573
VECTOR	pDONR223
NCBI SYMBOL	ASB17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtaaatctactaaattatgtgtaagactctgtccaagaagcaatatattctgcaatctctgacaaaattgttaaagaccctccctacagtttttgggtcagtgaggatcactgttacgaaccaaggattacagatcactggcaaaaattctgaggtatgtgacttggaattga
cgactactcacagattacattgcatgttggaataatcaggataccgttttgaagtaagtgttaacctcgacttcaatgaaatgtgtgaatacaattctgtactgggttttggcagaaaaggtaactcgtactgttggaattgcttcaagaagacaaaagactatgttcaagacagaagtgtga
acctggcactgatatggagaactttcacaccagtatactgccaagcccaatgaagtgccatcacacctctctttatgtagctcagacaagacagctcaatactcctcaaaaactactgcaatatggaatcttagaaaagagaaaaaacctatcaacattgtcttaacaatagactctacccctcga
gagtaagagtaattgtgatcgtgaattggctgacatccatgaagatgccaaaacatgtttgtactatgttccagagtcttctgtcatttcagtcaggaaataaagacacagctgagtttaggaagacatccaattatttcaaatgggttgattacattcctcaacaagatacaaatccatg
tgaactattacatcttgcagactaacatcaggatcaactattaaccaacaatatgtcccagatggaatatcttcaattctctcgtctacaaaactatctgaatttagaatactaa