

LOC400657-hypothetical LOC400657 Gene

Reference: PTXBC036588

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	long intergenic non-protein coding RNA 909		
For research use only.			

GENE DESCRIPTION	hypothetical LOC400657
GENE ID	400657
ORIGIN SPECIES	Human
ACCESSIONS	BC036588
VECTOR	pDONR223
NCBI SYMBOL	LOC400657
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcgcgtcgggtccaaggcaggcgtgtcagccgggaggggtgtctgggatttctagaggggtccagaacgggttctccaggccgacccgcctgtctctgcccggtctctccaccgcctccgcctctgttaaggaccccgagggcgccgcccgaacacggagaa
ggcccagaacggcgaggacccgagaggcctactgttgcaaacgcagcccgctgccctttgtctgcttccccggcgctggagtggtctctgtctgtctgtccgatttcgagaaatgactttcagctgtgcaaggcgagaatgctacttttagataa