

LEMD1-LEM domain containing 1 Gene

Reference: PTXBC036636

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT50; LEMP-1; LEM domain-containing protein 1; LEM domain protein 1; cancer/testis antigen 50		

For research use only.

GENE DESCRIPTION	LEM domain containing 1
GENE ID	93273
ORIGIN SPECIES	Human
ACCESSIONS	BC036636
VECTOR	pDONR223
NCBI SYMBOL	LEMD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggatgtgaagtgtctgactgtaattgcagaaccaacttgagaagcttggaatttcacctggcccaatactacctccaccagaaagtgtatgaaaaaagttagtacagttgttggtctcacctccctgtgcaccacctgtgatgaatggaccagagagctggatggagcgca
ggacagtgatgacagcgaagagcttaatacatctttgcaaggaaatcatcactctcaacaaaaaaagcaagaaactcaaaaatggcctgaggctccaccactaaacgcaaagctgtagataacctattgctggattataagcctccaagggaagaagggtggctgcaagagcaccaa
gcaccagaatcacatattgggactatcaccaagagagagactactgcgcggaagaccagactatcgagagctggagagaagaagggttccagtggtgctgaagcttgctgtgctgtatttcatcattgtggtgttctactgactgtggaataagtcgtgttgggttaa