



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

C19orf28-chromosome 19 open reading frame 28 Gene

Reference: PTXBC036706

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf28; PP3501; major facilitator superfamily domain-containing protein 12; major facilitator superfamily domain containing 12		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 28
GENE ID	126321
ORIGIN SPECIES	Human
ACCESSIONS	BC036706
VECTOR	pDONR223
NCBI SYMBOL	C19orf28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr