

CLDN12-claudin 12 Gene

Reference: PTXBC036754

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	claudin-12; claudin 12		
For research use only.			

GENE DESCRIPTION	claudin 12
GENE ID	9069
ORIGIN SPECIES	Human
ACCESSIONS	BC036754
VECTOR	pDONR223
NCBI SYMBOL	CLDN12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctgtcgggatgccacgcagccacagtcctttcctcctgtgtggaatcgctcagtagcaggcctcttgcaggagctctgttcccaactggagaaaattacgattgacacattcaacagaaacgagaagaacctgactgtttacacaggcctgtgggtgaaatgtgcccggtatg
acgggagcagtgactgctgatgtacgacactacttggtactcatcagttgaccagctggacctgcgtgtcctccagtttgcctaccctcagcagctgctgacgcatgggtgacctgtctctgctgattggaaatgtgcaacactgcttcaggtcctcgtgtcccaacatcaaatggc
caagtgtctgtcaatagtgaggttgcacactgtgtggtggtgctgtatcttctggcaggtactgtgagcctctcccatctatctgggtcatctttataacatccatctgaacaagaagttgagccagtcctttcatttgactatgcagtgatgtactatgtcctagtgctggggcctgttat
gacttccctatactatttattgtattgtacatgcaaatcttgccttctccttcttggcaaccattgtactccatccaccagtagtcatactactacagccctattcagcagctctcgcctctctgccattgaaattgacattccagtagtttcacacaccactaa