



Reference: PTXBC036755

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PDZRN1; ligand of Numb protein X 2; PDZ domain containing ring finger 1; PDZ domain-containing RING finger protein 1; numb-binding protein 2; ligand of numb-protein X 2		

For research use only.

GENE DESCRIPTION	ligand of numb-protein X 2
GENE ID	222484
ORIGIN SPECIES	Human
ACCESSIONS	BC036755
VECTOR	pDONR223
NCBI SYMBOL	LNK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]