



Reference: PTXBC036756

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G13; guanine nucleotide-binding protein subunit alpha-13; g alpha-13; guanine nucleotide binding protein (G protein), alpha 13; G protein subunit alpha 13		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), alpha 13
GENE ID	10672
ORIGIN SPECIES	Human
ACCESSIONS	BC036756
VECTOR	pDONR223
NCBI SYMBOL	GNA13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]